

Prodigy™

THROMBECTOMY SYSTEM

Product Specifications

Prodigy™ Catheters				
Product	Model	Length	Inner Diameter	Outer Diameter
Prodigy Catheter 5F	PRODIGY5F	160 cm	0.055 in	1.7 mm
Prodigy Catheter 6F	PRODIGY6F	137 cm	0.071 in	2 mm
Prodigy Catheter 8F	PRODIGY8F	117 cm	0.088 in	2.8 mm
Prodigy Catheter 8F-s	PRODIGY8F-S	50 cm	0.088 in	2.8 mm

Prodigy Twist™		
Product	Model	Length
Prodigy Twist 5F	TWIST5F	177 cm
Prodigy Twist 6F	TWIST6F	154 cm
Prodigy Twist 8F	TWIST8F	134 cm
Prodigy Twist 8F-s	TWIST8F-S	67 cm

TRUVIC™ Generator and Accessories	
Product	Model
TRUVIC Generator	PTG1B
TRUVIC Canister	PTC001B
TRUVIC Tubeset	PAT10261B

Hotshot™ Controller	
Product	Model
Hotshot Controller	HOTSHOT

Customer Service Information

Phone: +1-833-878-8422 | Fax: +1 650-618-2444
customerservice@imperativecare.com

