

# Product Specifications

Customer Service Information  
 CustomerService@ImperativeCare.com  
 +1 (408) 502-7548

## Symphony™

THROMBECTOMY SYSTEM

Symphony Catheters*	Model Number	Length (CM)	Sheath Compatibility
Symphony 16F Catheter, 82cm	SYMPHONY16F82	82	16F
Symphony 16F Catheter	SYMPHONY16F	117	16F
Symphony 24F Catheter	SYMPHONY24F	85	24F

\*All Symphony Catheters packaged with length-matched dilator

Symphony Auxiliary Devices	Model Number	Length (CM)
16F ProHelix™	PROHELIX16F	139
24F ProHelix	PROHELIX24F	102
Symphony Clot Container	SYMPHONYCC	N/A

Vascular Generator and Accessories	Model Number
Vascular Generator	PTG1B
Vascular Canister	PTCOO1B
Tubeset	PAT10261B



Rx only.

The Symphony 16F 82cm Thrombectomy System™ is intended for:

- The non-surgical removal of fresh, soft emboli and thrombi from blood vessels.
- Injection, infusion, and/or aspiration of contrast media and other fluids into or from a blood vessel.

The Symphony 16F 82cm Thrombectomy System is intended for use in the peripheral vasculature. It is not for use in the pulmonary vasculature.

The Symphony™ Thrombectomy System is intended for:

- The non-surgical removal of fresh, soft emboli and thrombi from blood vessels.
- Injection, infusion, and/or aspiration of contrast media and other fluids into or from a blood vessel.

The Symphony Thrombectomy System is intended for use in the peripheral vasculature. It is not for use in the pulmonary vasculature.

For complete product information, including indications, contraindications, warnings, precautions and adverse events, see product IFU included in product packaging,

visit [bit.ly/3pAaUlw](http://bit.ly/3pAaUlw), or contact Customer Service at:

**Phone:** +1-408-502-7548

**Fax:** +1-650-618-2444

**Email:** [customerservice@imperativecare.com](mailto:customerservice@imperativecare.com)

**Web:** [www.imperativecare.com/vascular](http://www.imperativecare.com/vascular)

Product image is not necessarily reflective of current branding.

© 2024 Imperative Care, Inc. All Rights Reserved. Imperative Care products are covered by U.S. and foreign patents. See [www.imperativecare.com/patents](http://www.imperativecare.com/patents) for a list of applicable U.S. patents. Unless otherwise specified, all product and service names appearing in this document are trademarks owned by or licensed to Imperative Care, Inc., its subsidiaries or affiliates. No use of any Imperative Care, Inc. trademark, trade name, or trade dress may be made without the prior written authorization of Imperative Care, Inc., except to identify the products or services of the company. See [www.imperativecare.com/trademarks](http://www.imperativecare.com/trademarks) for more information.

M10147.D