

Product Specifications


Customer Service Information
CustomerService@ImperativeCare.com
+1 (408) 502-7548

Symphony® Thrombectomy System

Symphony Catheter*	Model Number	Working Length (CM)	Overall Length (CM)	Sheath Compatibility (F)
16F Catheter, Short	SYMPHONY16F82	82	104	16
16F Catheter, Long	SYMPHONY16F	117	139	16
24F Catheter	SYMPHONY24F	85	103	24

*All Symphony Catheters packaged with length-matched dilator

Symphony Auxiliary Devices	Model Number	Length (CM)	Generator & Accessories	Model Number
24F Advance Long Dilator	ADVANCE	124	Generator	PTG1B
16F ProHelix	PROHELIX16F	139	Canister	PTC001B
24F ProHelix	PROHELIX24F	102	Tubaset	PAT10261B
Symphony Clot Containers	SYMPHONYCC	N/A		



Rx only.

The Symphony® Thrombectomy System is intended for:

- The non-surgical removal of fresh, soft emboli and thrombi from blood vessels.
- Injection, infusion, and/or aspiration of contrast media and other fluids into or from a blood vessel.

The Symphony Thrombectomy System is intended for use in the peripheral vasculature and for the treatment of pulmonary embolism.

The Symphony 16F 82cm Thrombectomy System shares the same peripheral indications described above; however it is not for use in the pulmonary vasculature.

For complete product information, including indications, contraindications, warnings, precautions and adverse events, see product IFU included in product packaging, visit bit.ly/3pAaUlw, or contact Customer Service at:

Phone: +1-408-502-7548

Fax: +1-650-618-2444

Email: customerservice@imperativecare.com

Web: www.imperativecare.com/vascular

Product image is not necessarily reflective of current branding.